

yuwell

POSITIVE AIRWAY PRESSURE UNITS

YH-830

- Specification
- Mode:CPAP、S、T、ST、VGPS
- Pressure:4-30cmH₂O
- Noise:28dBA±2dBA
- Efficient humidifier
- High Pressure Reminder
- Leak Reminder
- Asphyxiation Reminder
- Trigger Adjustment
- Cycle Adjustment
- Slope Adjustment
- Tidal Volume Adjustment
- Altitude Adjustment
- Leak Compensation
- Data Management
- Central Apnea Detection
- Smart Humidification



PACKAGING PARAMETERS

Model	QTN (PCS/CTN)	G.W/N.W (KG)	Size (cm)
YH-830	1	4.5/2.4	48.5×32.5×21.5

JIANGSU YUYUE MEDICAL EQUIPMENT & SUPPLY CO., LTD

Manufacturing Address: No.1 Baisheng Road, Development Zone, Danyang, Jiangsu Province, P.R.China

Tel: 86-0511-8690 0859 86-0511-8690 0860

E-mail: yuwellsales@yuyue.com.cn; sales@yuyue.com.cn

www.yuwell.com

f | t | in | @ | Yuwell Global